

Analysis of sodium bicarbonate present in a commercial sample of soda mint tablet.

By Dr. Chitra Joshi

Aim: To determine the amount of sodium bicarbonate present in the given sample of soda mint tablet.

Theory:

Soda mint tablet contains sodium bicarbonate which is a weak base and, therefore, soda mint tablets are used as antacids. The amount of sodium bicarbonate present in the tablet can be estimated by volumetric analysis. This is an acid-base titration. Hydrochloric acid neutralizes Sodium bicarbonate as follows,



Methyl orange is used as an indicator in this titration to determine the neutralization point. NaHCO_3 present in a given sample tablet can be calculated using the relation,

1 M 1000 cm^3 HCl \equiv 84 g anhydrous NaHCO_3

Volumetric analysis

Volumetric analysis is a technique that employs the measurement of volumes to quantitatively determine the amount of a substance in solution. In any reaction between two or more species, the reaction equation shows the stoichiometric ratio of the reacting species. Hence, if the concentration of one of the solutions is known, the concentrations of the others can be determined from the volumes used.

Succinic acid ($\text{C}_4\text{H}_6\text{O}_4$, Molar mass 118g/mole) is used to standardize alkali (NaOH).

Succinic acid is a stable solid organic acid, known molar solution of it can be prepared very easily. Using standardized alkali (NaOH) other acid (say HCl) can be standardize. Since the relation 1 M 1000 cm^3 HCl \equiv 84 g anhydrous NaHCO_3 is given, standardization of HCl is to be done carefully and accurately.

Requirements:

No.	Name of the apparatus/Chemical	Quantity
1	HCl solution	100 cm^3
2	0.1M(aprox..)NaOH	100 cm^3
3	Soda mint tablets.	3
4	Phenolphthalein indicator.	1 small indicator bottle.
5	Methyl orange indicator.	1 small indicator bottle.
6	Distilled water (DW)	1 bottle
7	Burettes (25 cm^3)	2 (B and A)
8	Syringes (10 cm^3)	3 (N, H and S)
9	Volumetric flasks(100 cm^3)	1 (S)
10	Beaker (100 cm^3)	1 (S)
11	Conical flasks	3
12	Stirrer	1
13	Dropper	1
14	Funnels	2
15	0.05M Succinic acid solution	100 cm^3

Procedure:**A. Preparation of Soda mint solution:**

Place the supplied three soda mint tablets in the beaker 'N'. Add approximately 20 cm³ distilled water to it crush and stir well to dissolve it completely (Solution will be slightly turbid due to undissolved binders). Filter this solution through cotton wool (**Plug the funnel with cotton to remove undissolved binders**). Collect the filtrate and washings of beaker 'N' in 100 cm³ beaker 'S'. Transfer the filtrate quantitatively to 100cm³ standard volumetric flask 'S', then add distilled water to the volumetric flask till mark on it (if required use dropper) to make the volume of the solution 100 cm³. Shake the solution to homogenize.

B. Standardization of 0.1M (approx.) NaOH solution.

1. Wash the burette (B) thoroughly with distilled water. Rinse and fill it with the supplied 0.1M (approx.) NaOH solution.
2. Using syringe (N) take out 10 cm³ of 0.05M succinic acid solution in a 100cm³ conical flask. Add two drops of phenolphthalein indicator to it.
3. Titrate the contents of the flask against the supplied NaOH solution from the burette, with constant swirling till light pink colour is obtained.
4. Repeat the procedure to get constant readings.
5. Enter your readings in the observation table I

Observations Table I:**Titration of 0.05 M Succinic acid versus NaOH solution**

Readings	Pilot	I	II	III	CBR
Final	----- to -----				B cm ³
Initial		0.0	0.0	0.0	
Difference					

Calculation:

Standardization of NaOH solution

$$M_{\text{NaOH}} \times V_{\text{NaOH}} = M_{\text{Succinic acid}} \times V_{\text{Succinic acid}} \times 2 \quad (\text{since succinic acid is dibasic acid})$$

$$M_{\text{NaOH}} = (0.05 \times 10 \times 2) \div \mathbf{B}$$

$$M_{\text{NaOH}} = \mathbf{X} \text{ M}$$

Exact molarity of NaOH solution is **X M**

C. Determination of exact molarity of the supplied HCl.

1. Fill the burette (B) with the **X M** NaOH solution.

- Using the syringe (H) take out 10 cm³ of the supplied HCl acid solution in a 100cm³ conical flask. Add two drops of phenolphthalein indicator to it.
- Titrate the contents of the flask against the X M NaOH solution from the burette, with constant swirling till light pink colour is obtained.
- Repeat the procedure to get constant readings.
- Enter your readings in the observation table II

Observations Table II:

Titration of X M NaOH solution versus HCl solution

Readings	Pilot	I	II	III	CBR
Final	-----				C cm ³
Initial	to	0.0	0.0	0.0	
Difference	-----				

Calculation:

Standardization of HCl solution

$$M_{\text{HCl}} \times V_{\text{HCl}} = M_{\text{NaOH}} \times V_{\text{NaOH}}$$

$$M_{\text{HCl}} = (X \times C) \div 10$$

$$M_{\text{HCl}} = Y M$$

Exact molarity of the HCl solution is Y M

D. Titration of soda mint solution with standard Y M HCl solution.

- Wash the burette (A) thoroughly with distilled water. Rinse and fill this burette (A) with standard Y M HCl solution.
- Transfer the soda mint tablet solution in the container 'S'.
- Using syringe (S) take out 10cm³ of the soda mint solution in 100cm³ conical flask. Add four drops of methyl orange indicator to it.
- Titrate the contents of the flask against standard Y M HCl solution from the burette, with constant swirling till orange colour is obtained.
- Repeat the procedure to get constant readings.
- Enter your readings in the observation table III

Observations Table III:**Titration of soda mint solution versus Y M HCl solution**

Readings	Pilot	I	II	III	CBR
Final	-----				D cm³
Initial	to	0.0	0.0	0.0	
Difference	-----				

Calculation:

Amount of Sodium bicarbonate present in 10 cm³ diluted solution of soda mint tablet:

Using the relation

$$1 \text{ M} \quad 1000 \text{ cm}^3 \text{HCl} \quad \equiv \quad 84 \text{ g anhydrous NaHCO}_3$$

$$Y \text{ M} \quad D \text{ cm}^3 \text{HCl} \quad \equiv \quad (Y \times D \times 84) \div 1000 = Z \text{ g anhydrous NaHCO}_3$$

100 cm³ diluted solution of soda mint tablets contain **10 Z** g anhydrous NaHCO₃

Since three tablets dissolved in 100cm³ solution, **10Z** g of anhydrous NaHCO₃ is present in 3 tablets hence one soda mint tablet contain average (**10Z** ÷ 3) g anhydrous NaHCO₃

1) Exact molarity of HCl

2) Amount of anhydrous NaHCO₃ present in the given sample solution

g

3) Average amount of anhydrous NaHCO₃ in the given sample of one soda mint tablet

g

Observations:

I) Reading of part B (NaOH \propto Succinic acid)=12.5 cm³

II) Reading of part C (NaOH \propto HCl acid)=12.4 cm³

III) Reading of part D (HCl \propto NaHCO₃)=10.5 cm³

Calculations:

I) Exact molarity of NaOH is 0.08M

II) Exact molarity of HCl is 0.0992M = 0.1M

III)

1 M 1000 cm³ HCl \equiv 84 g anhydrous NaHCO₃

0.1M 10.5 cm³ HCl = (0.1 x 10.5 x 84) / 1000 = 0.882g

Therefore for 100cm³ of soda mint solution amt. of NaHCO₃ = 0.882g

For three tablets amt. of NaHCO₃ = 0.882g

For one tablet amt. of NaHCO₃ = 0.882/3 = 0.294g

In the commercial sample of soda mint tablet amount of NaHCO₃ is 0.300g
